





Advantages





- Real security of results with the COBAS INTEGRA clot detection system
- cobas c packs make reagent preparation a task of the past
- Workflow consolidation
 with the choice of over 140 applications
- Sophisticated software programs to simplify operation and training
- Reduced costs through increased efficiency
- Continuous access
 to up to 90 samples for a smoother workflow

COBAS INTEGRA 400 plus is designed to consolidate testing and increase efficiency while reducing total running costs of the laboratory.

The COBAS INTEGRA 400 plus system is the right choice for consolidation in the low workload laboratory and special chemistry testing in medium volume sites.

The COBAS INTEGRA systems' extensive assay menu and innovation of the **cobas c** pack are combined with four proven measurement technologies and sophisticated, easy-to-use software.



Convenience



cobas c pack

Patented reagent carrier design plus on-board refrigeration prevents reagent evaporation and degradation, ensuring long-term on-board stability and long calibration intervals.

Each **cobas c** pack holds all necessary reagents for up to 800 determinations, maximising on-board capacity and reducing the need for operator intervention and frequent reloading.

Reagent preparation is a task of the past. The barcoded **cobas c** packs are automatically handled by the system, reducing the possibility of errors and saving staff time.

Compact design of the reagent carrier means that more than a year's supply of tests for a typical laboratory can easily be stored in a single refrigerator.

Continuous sample access

Sample racks allow continuous access as tests are completed and new samples arrive, eliminating unnecessary test delays.

Up to 90 barcoded samples, in six 15 position racks, can be loaded on-board. STATs are automatically given priority and immediately processed once the barcode is read to speed up delivery of results to the physician.

Serum, plasma, urine, hemolysate, whole blood and CSF can be used in barcoded primary or secondary tubes. Microcups can be used for low volume samples, e.g. paediatric and CSF.

Workflow consolidation

More than 140 applications for all types of sample matrices are measured with one of 4 different measuring technologies – absorbance photometry, turbidimetry, fluorescence polarisation and ion selective potentiometry.


The COBAS INTEGRA 400 plus system has up to 36 different on-board assays. For quick access to additional tests, **cobas c** packs can be exchanged while the system is running. Off-board storage of pre-calibrated **cobas c** packs allows full flexibility.

The broad menu includes clinical chemistry, specific proteins, therapeutic drug monitoring, drug of abuse testing and much more. Consolidation on a single system improves turnaround times by reducing the need for sample splitting.

On-board storage of 32 **cobas c** packs at 12°C ensures reagent stability up to 3 months, reducing reagent waste and costs.

For high volume tests, multiple ${\bf cobas}\ {\bf c}$ packs of the same assay can be on-board for automatic switch-over during run.





Reliable and efficient

Sophisticated software

On-screen inventory management tracks reagents available on board and forecasts if additional reagent packs are required based on daily workload statistics.

The intelligent processing on COBAS INTEGRA 400 plus automatically schedules each sample's tests to minimise analyser turnaround time and increase throughput, while allowing immediate STAT testing 24 hours a day.

Quality control programs monitor precision and accuracy to ensure reliable results.

Validated sample results are continuously transmitted to the LIS via the host query communication mode, guaranteeing fast result availability.

The software automatically tracks service needs and alerts the user when maintenance is due. Activities are recorded in the electronic log to support compliance with Good Laboratory Practice and Accreditation.

Remote diagnostics by modem allows Roche service personnel to resolve errors quickly and efficiently without an on-site visit.



Reduced costs through increased efficiency

Consolidation eliminates pre-analytical phase tasks by reducing sample splitting to improve workflow and decrease handling costs.

Long on-board reagent stability and fewer calibrations mean less waste and lower cost per reportable result.

A single workstation optimises use of time and budget by reducing costs of operating, maintaining and servicing multiple analytical systems.

Processing time can be reduced by up to 50% compared with dedicated analysers.


Post-analytical tasks such as sample tracking, rerunning, archiving and disposal are greatly simplified by use of a single consolidated system – COBAS INTEGRA 400 plus.

Real security

Highly sensitive pressure sensors detect incorrect pipetting, even at 2 µL sample volume.

Handling of clotted samples is improved with immediate flagging of faulty samples.

System wash pressures have been optimised to fully remove the clot from the sample probes, avoiding any maintenance actions and interruptions to workflow.



System overview



Automatic start-up

- Configurable start-up options ensure that the system is ready whenever needed
- Access to the full COBAS INTEGRA reagent menu 24 hours a day guarantees quick after hours turn around time
- Automatic system maintenance significantly reduces operator interventions

2 Result security

- Highly sensitive pressure sensors detect incorrect pipetting, even at 2 µL sample volume
- Handling of clotted samples is improved with immediate flagging of faulty samples
- System wash pressures have been optimised to fully remove the clot from the sample probe, avoiding maintenance actions and interruptions to workflow

3 ISE modul

• Sodium, Potassium, Chloride and Lithium

Flexible change of onboard tests

- Onboard capacity of 8 racks with 4 cobas c packs each
- Easy change of menu via the reagent racks
- Refrigerated storage for long calibration stability





5 Continuous sample loading

- 6 sample racks for various primary and secondary cups
- Up to 90 samples onboard, racks can be replaced after sample pipetting
- Positive sample identification via barcode

6 Analytical unit


- Precise temperature controlled reaction rotor ensure stable test conditions and reliable results
- Sensitive photometric and turbidimetric measurements
- Extremely stable LED light source for fluorescence polarisation

2 Enhanced user interface

- System software is object-oriented, based on Windows XP, for standardised ease-ofuse
- Intuitive online help provides instant and precise multilingual support
- Easy to read 19 inch flatscreen monitor

Specifications

System	Random and continuous access, samp Integration of 4 measuring principles 36 tests onboard Absorbance photometry: Turbidimetry: Fluorescence polarimetry: Ion-selective electrode potentiometry:	Enzymes and substrates Specific proteins, drugs of abuse Therapeutic drugs, thyroid tests
Test throughput	Up to 400 tests/hr (including ISE)	
Sample types	Serum, plasma, urine, CSF, hemolysate and whole blood (HbA1c)	
Sample handling	90 primary or secondary tubes on-board Up to 6 x RD15 racks on-board Cooled rack position for controls and calibrators Automatic sample dilution and concentration Barcode reading via laser scanner, with immediate STAT recognition	
Sample container types	Primary tubes: 5 to 10 mL; 16x100, 16x75, 13x100, 13x75 mm with closed cap for Hitachi microcup, 500 μL; Hitachi standard cup, 1.5 mL; cobas cup, 650 μL; Eppendorf cup, 1.5 mL; Cup on tube definable	
Sample volume	Typically 2 to 10 μL per test, ISE indirect 20 μL, ISE direct 97 μL	
Sample barcode types	Code 128, Codabar, Code 2 of 5 interleaved, Code 39	
Reagents onboard capacity	32 cobas c packs, 50-800 tests per cobas c pack Up to 8 racks of 4 cobas c packs on-board Reagent compartment cooled to 10-15°C Onboard stability up to 3 months, calibration typically each lot	
Reaction cells	Holds 1,000 disposable cuvettes with 5 mm path length, and 120-240 µL reaction volume	
Control unit	HP workstation running Windows XP Intel Pentium IV with 512 MB RAM Dual 40 GB hard drives, CD ROM, floppy drive and inbuilt modem	
System interfaces	RS 232 serial interface, bi-directional, inbuilt modem for remote diagnostics access	
Electrical requirements	100-125/200-240 Volts AC, 50 or 60 Hz, consumption 1,200 VA	
Physical dimensions	Width: 135 cm (53.1 in) Depth: 66 cm (25.9 in) Height: 75 cm (28.5 in)	
Weight	320 kg (506 lbs)	
Water requirements	Up to 2 liters per hour in operating mode, Type 1 NCCLS	
Certification	CE, UL, C-UL	



COBAS, COBAS C, COBAS INTEGRA and LIFE NEEDS ANSWERS are trademarks of Roche. PENTIUM is a trademark of Intel. WINDOWS is a trademark of Microsoft Corporation.

©2009 Roche

Roche Diagnostics Ltd. CH-6343 Rotkreuz Switzerland www.roche.com

Enzymes 5'Nucleotidase²

Test menu

Substrates
Albumin
Albumin BCP ²
Ammonia
Bicarbonate
Bilirubin - direct
Bilirubin - total
Calcium
Cholesterol
HDL Cholesterol Direct
LDL Cholesterol Direct
Creatinine
Fructosamine
Glucose
Homocysteine ²
Iron
Lactate
Magnesium
Phosphorus
Total Protein
Protein Urine / CSF
Triglycerides
Total Bile Acids ²
UIBC
Urea / BUN
Uric Acid
Zinc ²

ACE ²
Acid Phosphatase
Alkaline Phosphatase
ALT / GPT
Amylase - pancreatic
Amylase - total
Antistreptodornase ²
AST / GOT
Cholinesterase / Dibucain
CK
CK-MB
GGT
GLDH
HBDH
_DH (L-P)(P-L)
_ipase
SE
Chloride
_ithium
Potassium
Sodium

Specific proteins
α1-Acid Glycoprotein
α1-Antitrypsin
α1-Microglobulin (c pack MULTI)
α2-Macroglobulin²
β2-Microglobulin²
Albumin (Turbidimetric)
APO A1
АРО В
ASLO
C3c
C4
Ceruloplasmin
CRP
Cystatin C
Ferritin
Haptoglobin
HbA1c (whole Blood, Hemolysate)
hsCRP
lgA
lgG
lgM
IMA²
Kappa (c pack MULTI)
Lambda (c pack MULTI)
Lp (a)
Myoglobin
Prealbumin
Rheumatoid Factor
Soluble Transferrin Receptor
Transferrin





Drugs of abuse Amphetamines Barbiturates - urine/serum Benzodiazepines - urine/serum Buprenorphine² Cannabinoids Cocaine Metabolite EDDP¹ Ethanol - urine/serum Lysergic Acid Diethylamide Methadone Methaqualone Opiates Oxycodone² Phencyclidine Propoxyphene Tricyclic Antidepressants²

Therapeutic drug monitoring		
Acetaminophen		
Amikacin		
Carbamazepine		
Caffeine ²		
Cyclosporine		
Digitoxin		
Digoxin		
Ethosuximide ²		
Gentamicin		
Lithium		
Lidocaine		
Methotrexate ²		
MPA (total)		
NAPA		
Netilmycin ²		
Phenobarbital		
Phenytoin – free/total		
Primidone		
Procainamide		
Quinidine		
Salicylate		
Sirolimus ²		
Tacrolimus ²		
Theophylline		
Tobramycin		
Valproic Acid -free/total		
Vancomycin		

Others	
AT III	
D-Dimer	
Factor XIII ²	
Serum Index	
T4	
T-Uptake	

Note:

Not all tests may be available in all countries due to local registration requirements. Please contact your local nearest Roche office for full details.

Urine and CSF applications are available for selected tests.

- ¹ Assays currently under development
- ² Available as Development Channel

Last review April 2009

COBAS, COBAS INTEGRA and LIFE NEEDS ANSWERS are trademarks of Roche.

©2009 Roche

Roche Diagnostics Ltd. CH-6343 Rotkreuz Switzerland www.roche.com